




# Test Certificate

Part 1 Liter Vessel inactinic brown  
 Part Number 02503-2  
 Serial Number 84112 - 87537



		USP / Ph.Eur. Tolerances	
<b>Dimensions</b>	Vessel	Inside Ø a	100 mm 98 – 106 mm
		Height b	200 mm 160 – 210 mm
		With hemispherical bottom	
<b>Material</b>	Borosilicate glass inactinic brown		

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopoeia requirements (at date of approval, see below).

Aesch, 2021-11-02

Bruno Balli  
 SOTAX Quality Management  
 SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland