




Part 1 Liter Vessel  
 Part Number 14143-11  
 Serial Number P35 001172 – P35 001274



# Test Certificate

		USP / Ph.Eur. Tolerances	
<b>Dimensions</b>	Vessel	InsideØ a	102.9 to 105.1 mm
		Height b	168.2 mm
		With hemispherical bottom	160 – 210 mm
<b>Material</b>	Borosilicate glass		

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopoeia requirements (at date of approval, see below).

Navi Mumbai, 2021-12-22

12/22/2021 2:08 PM

Moreshwar Pulekar  
 SOTAX Quality Management  
 C-64 | First Floor | TTC Industrial area | Turbhe | Navi Mumbai 400705 | India