




Part 1 Liter Vessel  
 Part Number 14143-11  
 Serial Number P35 001846 – P35 001958



		USP / Ph.Eur. Tolerances		
<b>Dimensions</b>	Vessel	Inside Ø a	102.9 to 105.1 mm	98 – 106 mm
		Height b	200 mm	160 – 210 mm
		With hemispherical bottom		
<b>Material</b>	Borosilicate glass			

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Navi Mumbai, 2022-12-06

06-12-2022 15:49:26

Moreshwar Pulekar  
 SOTAX Quality Management  
 C-64 | First Floor | TTC Industrial area | Turbhe | Navi Mumbai 400705 | India

# Test Certificate