


Test Certificate

Part: **1 Liter Vessel**
 Number: 14322-01
 Serial Number: 73994 - 74255

For SOTAX AT Dissolution Testing Systems



Dimensions:	Inside Ø	a	98.5 to 101.5 mm	USP/Ph.Eur. Tolerances	98 – 106 mm
	Height	b	200 mm		160 – 210 mm
	With hemispherical bottom				
Material:	Borosilicate glass				

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at time of delivery).

SOTAX Quality Management
 Bruno Balli
 2018-03-14


 Signature