




Part China Vessel 250 mL
 Part Number 14536-01
 Serial Number 79206 - 80255



Test Certificate

			PPRC Tolerances	
Dimensions	Vessel	Inside Ø a	62 mm	± 3 mm
		Height b	126 mm	± 6 mm
		With hemispherical bottom		
Material	Borosilicate glass			

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with PPRC Pharmacopeia requirements (at date of approval, see below).

Aesch, 2019-10-03

Bruno Balli
 SOTAX Quality Management
 SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland