




Part 1 Liter Vessel
 Part Number 14537-04
 83672Serial Number 83536 - 83541



			Tolerances	
Dimensions	Vessel	Inside Ø a	102.9 to 105.1 mm	--
		Height b	155.1 mm	± 0.7 mm
		With hemispherical bottom		
Material	Borosilicate glass			

Aesch, 2020-09-14

Bruno Balli
 SOTAX Quality Management
 SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland

Test Certificate