

Test Certificate

Part: **Nozzle 10 mm**
Number: 231202

For SOTAX PF1



			USP/Ph.Eur. Tolerances
Dimensions:	Upper Inside Ø	30.2 mm	0 / - 0.1 mm
	Bottom Inside Ø	10.00 mm	± 0.01 mm
	Chamfer	1.0 mm	± 0.2 mm
	Funnel (Cone)	40 °	± 0.8 °

Material: **Stainless Steel (electrolytic polished)**

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at time of delivery).

SOTAX Quality Management
Bruno Balli
2018-08-13


Signature