




Part Plastic Sieve 10 mesh (6 Tubes)  
 Part Number 23440-01



**USP / Ph.Eur.  
Tolerances**

<b>Dimensions</b>	A	Aperture size	1.8 to 2.2 mm	--
		Aperture shape	0.57 to 0.66 mm	--
<b>Material</b>	Plastic RTP 702 - HDPE			

Remark: Aperture shape square

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP 701 and Ph. Eur. 2.9.1 Pharmacopeia requirements (at date of approval, see below).

Aesch, 2021-02-25

Bruno Balli  
 SOTAX Quality Management  
 SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland

# Test Certificate